

Important information for the handling of the urine collection vessel

Avoid any spilling by following this instruction:

- Remove the spout cover from the green lid
- Pull it in direction of the arrows
▶▶ (outside in)



- Keep bottle and vessel at an angle
- The spout of the collection vessel must rest on the rim of the bottle.



- Avoid any distance



 **berlingerspecial**
feel safe

Phone: +41 71 982 88 11
E-mail: info.specialberlinger.ch

The latest versions of the user manual and the instructional video can be found on our website: www.berlinger.ch